



Product Datasheet

CD147 / Emmprin / Basigin Antibody (orb1824859)


Description	BSG is expressed more intensely on thymocytes than on mature peripheral blood T cells. It is importa...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Recombinant full-length human protein was used as the immunogen for the CD147 antibody.
Storage	Aliquot the CD147 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	P35613
Hazard Information	This CD147 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon carcino...



IHC staining of FFPE human placental tis...



IHC staining of FFPE human renal cell ca...