



Product Datasheet

NeuN Antibody / Fox3 / Rbfox3 (orb1824937)


Description	NeuN antibody specifically recognizes the DNA-binding, neuron-specific protein NeuN, which is present...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the RBFOX3 antibody.
Storage	Aliquot the RBFOX3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	A6NFN3
Hazard Information	This RBFOX3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human brain tissue ...



IHC staining of FFPE human cerebrum with...



SDS-PAGE analysis of purified, BSA-free ...