



Product Datasheet

NeuN Antibody / Fox3 / Rbfox3 (orb1824939)


Description	This antibody specifically recognizes the DNA-binding, neuron-specific protein NeuN/Rbfox3, which is...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the RBFOX3 antibody.
Storage	Aliquot the RBFOX3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	A6NFN3
Hazard Information	This RBFOX3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human brain tissue ...



IHC staining of FFPE human placental tis...



IHC staining of FFPE human stomach tissu...