



Product Datasheet

CD31 Antibody / PECAM-1 (orb1824983)

Description	CD31 (PECAM-1) is a transmembrane glycoprotein member of the immunoglobulin supergene family of adhe...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 538-738) from the human protein was used as the immunogen for the CD31 antibody.
Storage	Aliquot the CD31 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG. kappa
Clonality	Recombinant
Uniprot ID	P16284
Hazard Information	This CD31 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human uterus tissue...



IHC staining of FFPE human placental tis...