



Product Datasheet

MPO Antibody / Myeloperoxidase (orb1825008)

Description	Myeloperoxidase (MPO) also called the peroxidase (POD), is an important marker of bone marrow cells....
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human MPO protein was used as the immunogen for the MPO antibody.
Storage	Aliquot the MPO antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P05164
Hazard Information	This MPO antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human spleen tissue...



SDS-PAGE analysis of purified, BSA-free ...