

---

## Product Datasheet

**NME2 Antibody / Nucleoside diphosphate kinase B  
/ NDKB (orb1825028)**

**Description**

The nm23 gene, a potential suppressor of

hybri...

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P, WB

**Immunogen**

Recombinant full-length human protein was used as the immunogen for the NDKB antibody.

**Storage**

Aliquot the NDKB antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Formula**

1 mg/ml in 1X PBS; BSA free, sodium azide free

**Isotype**

Mouse IgG1, kappa

**Clonality**

Monoclonal

**Uniprot ID**

**P22392**

**Hazard Information**

This NME2 antibody is available for research use only.

**Dilution Range**


Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date**


12 months from date of receipt.




Western blot testing of human MCF-7 cell...



IHC staining of FFPE human kidney tissue...



IHC staining of FFPE human kidney tissue...