

Product Datasheet

**NME2 Antibody / Nucleoside diphosphate kinase B
/ NDKB (orb1825029)**

Description

The nm23 gene, a potential suppressor of

hybri...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant full-length human protein was used as the immunogen for the NDKB antibody.

Storage

Aliquot the NDKB antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Uniprot ID

P22392

Hazard

Information

This NME2 antibody is available for research use only.

Dilution Range


Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



Western blot testing of human MCF-7 cell...



IHC staining of FFPE human kidney tissue...



IHC staining of FFPE human kidney tissue...