

Product Datasheet

NME2 Antibody / Nucleoside diphosphate kinase B / NDKB (orb1825034)

Description

The nm23 gene, a potential suppressor of

hybri...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant full-length human NME2 protein was used as the immunogen for the NME2 antibody.

Storage

Aliquot the NME2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Mouse IgG2, kappa

Clonality

Monoclonal

Uniprot ID

[P22392](#)

Hazard

Information


This NME2 antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



IHC staining of FFPE human skin tissue w...



Western blot testing of human MCF-7 cell...



Western blot testing of human HEK293 cel...