



Product Datasheet

NME2 Antibody / Nucleoside diphosphate kinase B / NDKB (orb1825035)


Description	The nm23 gene, a potential suppressor of metastasis, was originally identified by differential hybri...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human NME2 protein was used as the immunogen for the NME2 antibody.
Storage	Aliquot the NME2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	P22392
Hazard Information	This NME2 antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...



IHC staining of FFPE human skin tissue w...



Western blot testing of human MCF-7 cell...