

Product Datasheet

**NME2 Antibody / Nucleoside diphosphate kinase B
/ NDKB (orb1825036)**

Description

The nm23 gene, a potential suppressor of

hybri...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

Applications

IHC-P, WB

Immunogen

Recombinant full-length human NME2 protein was used as the immunogen for the NME2 antibody.

Storage

Aliquot the NME2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2, kappa

Clonality

Monoclonal

Uniprot ID

P22392

Hazard

Information


This NME2 antibody is available for research use only.

Dilution Range


Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



Western blot testing of human HEK293 cel...



IHC staining of FFPE human skin tissue w...



Western blot testing of human MCF-7 cell...