
Product Datasheet

GFAP Antibody / Glial Fibrillary Acidic Protein (orb1825065)

Description

This mAb recognizes a protein of ~50kDa which is

Species/Host

Mouse

Reactivity

Bovine, Canine, Equine, Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

A recombinant partial protein sequence (within amino acids 150-250) from the human protein was used as the immunogen for the GFAP antibody.

Storage

Aliquot the GFAP antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Uniprot ID

P14136

Hazard Information


This GFAP antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date


12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



IHC staining of FFPE human brain tissue ...



IHC staining of FFPE dog brain with GFAP...