

Product Datasheet

C5R1 (phospho-Ser338) antibody (orb182512)

Description	C5R1 (phospho-Ser338) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human C5R1 around the phosphorylation site of Ser338 (TR(p-S)TV)
Target	phospho-C5R1 (Ser338)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	P21730
Expiration Date	12 months from date of receipt.