



Product Datasheet

p53 Antibody / TP53 (orb1825265)

Description	Recognizes a 53kDa protein, which is identified as p53 suppressor gene product. It reacts with the m...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human full-length protein was used as the immunogen for the TP53 antibody.
Storage	Aliquot the TP53 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P04637
Hazard Information	This TP53 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder tissue...



SDS-PAGE analysis of purified, BSA-free ...