



Product Datasheet

MMP2 Antibody (orb1825290)

Description	It recognizes a protein of 72kDa, which is identified as MMP2. The matrix metalloproteinases (MMP) a...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human protein was used as the immunogen for the Matrix metalloproteinase 2 antibody.
Storage	Aliquot the Matrix metalloproteinase 2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P08253
Hazard Information	This Matrix metalloproteinase 2 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human placental tis...



SDS-PAGE analysis of purified, BSA-free ...