



## Product Datasheet

# TPSAB1 Antibody / Mast Cell Tryptase (orb1825298)

<b>Description</b>	Tryptases comprise a family of trypsin-like serine proteases (peptidase family S1). Tryptases are st...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	Recombinant human Tryptase 1 protein fragment (within amino acids 100-250) was used as the immunogen for the recombinant TPSAB1 antibody.
<b>Storage</b>	Aliquot the recombinant TPSAB1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	1 mg/ml in 1X PBS; BSA free, sodium azide free
<b>Isotype</b>	Mouse IgG1, kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q15661</a>
<b>Hazard Information</b>	This recombinant TPSAB1 antibody is available for research use only.
<b>Dilution Range</b>	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
<b>Expiration Date</b>	12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



SDS-PAGE analysis of purified, BSA-free ...