



Product Datasheet

Major Vault Protein Antibody / MVP (orb1825367)

Description	Recognizes a protein of 104kDa-110kDa, characterized as major vault protein (MVP). Vaults are large ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Proteins precipitated from human breast cancer MCF-7 cells were used as the immunogen for the recombinant MVP antibody.
Storage	Aliquot the recombinant MVP antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	Q14764
Hazard Information	This recombinant MVP antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate tiss...



Western blot testing of human U-87, A549...