
Product Datasheet

**CD19 Antibody / B-lymphocyte marker
(orb1825370)**

Description

CD19 is a transmembrane glycoprotein that contains


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

IHC-P



IHC staining of FFPE human tonsil tissue...

Immunogen

A recombinant fragment of human CD19 protein (around aa1-300) was used as the immunogen for the CD19 antibody.

Storage

Aliquot the CD19 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Mouse IgG1, kappa

Clonality


Monoclonal

Uniprot ID

[P15391](#)

Hazard Information

This CD19 antibody is available for research use only.




IHC staining of FFPE human tonsil tissue...

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...