
Product Datasheet

**CD19 Antibody / B-lymphocyte marker
(orb1825371)**

Description

CD19 is a transmembrane glycoprotein that contains


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

IHC-P



IHC staining of FFPE human tonsil tissue...

Immunogen

A recombinant fragment of human CD19 protein (around aa1-300) was used as the immunogen for the CD19 antibody.

Storage


Aliquot the CD19 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide



IHC staining of FFPE human tonsil tissue...

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Uniprot ID


P15391

Hazard Information

This CD19 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT



Analysis of a HuProt (TM) microarray con...

Expiration Date

12 months from date of receipt.