

Product Datasheet

CENTB1 (phospho-Ser554) antibody (orb182538)

Description	CENTB1 (phospho-Ser554) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CENTB1 around the phosphorylation site of Ser554 (PG(p-S)LR)
Target	phospho-CENTB1 (Ser554)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81 kDa
Uniprot ID	Q15027
Expiration Date	12 months from date of receipt.