



Product Datasheet

GCHFR Antibody / GFRP (orb1825430)

Description	GTP cyclohydrolase I feedback regulatory protein (GFRP) is encoded by the gene GCHFR. GFRP mediates ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human protein was used as the immunogen for the GCHFR antibody.
Storage	Aliquot the GFRP antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	P30047
Hazard Information	This GFRP antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human liver in HCC ...



SDS-PAGE analysis of purified, BSA-free ...