

Product Datasheet

IgA Heavy Chain Antibody (orb1825439)

Description	This mAb is specific to heavy chain of IgA and shows minimal cross-reaction with heavy chains of oth...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Full-length native purified protein corresponding to Human IgA was used as the immunogen for the recombinant IgA antibody. It is reactive with both IgA1 and IgA2 subclasses of Alpha heavy chain.
Storage	Aliquot the recombinant IgA antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P01877, P01876
Hazard Information	This recombinant IgA antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.