



Product Datasheet

SNW1 Antibody / NCoA-62 / SKIP (orb1825525)


Description	Nuclear receptor coactivator NCOA62/SNW1 (also called Nuclear Protein SkiP, SKIIP, Ski-interacting p...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF
Immunogen	A recombinant fragment of human SNW1 protein was used as the immunogen for the SNW1 antibody.
Storage	Aliquot the SNW1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a
Clonality	Monoclonal
Uniprot ID	Q13573
Hazard Information	This SNW1 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 2-4ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescent staining of PFA-fixed ...



Flow cytometry testing of PFA-fixed huma...



SDS-PAGE analysis of purified, BSA-free ...