



Product Datasheet

NKX2.1 Antibody / TTF-1 (orb1825555)


Description	Recognizes a protein of 40kDa, identified as Thyroid transcription factor-1 (TTF-1). TTF-1 is a memb...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human TTF-1 protein was used as the immunogen for the TTF-1 antibody.
Storage	Aliquot the TTF-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P43699
Hazard Information	This TTF-1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thyroid tissue...



IHC staining of FFPE human thyroid tissue...



SDS-PAGE analysis of purified, BSA-free ...