

Product Datasheet

IgD Heavy Chain Antibody (orb1825564)

Description Immunoglobulins are four-chain, Y-shaped, monomeric structures comprised of two identical heavy chain...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC-P

Immunogen A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the recombinant IgD Heavy Chain antibody.

Storage Aliquot the recombinant IgD Heavy Chain antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note For research use only

Formula 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype Rabbit IgG, kappa


Clonality Recombinant

Uniprot ID **P01880**


Hazard Information This recombinant IgD Heavy Chain antibody is available for research use only.

Dilution Range Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date 12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...