


Product Datasheet

MYH7 Antibody / Myosin 7 (orb1825579)

Description	Myosin heavy chains are ubiquitous Actin-based motor proteins that convert the chemical energy deriv...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human protein was used as the immunogen for the Myosin 7 antibody.
Storage	Aliquot the Myosin 7 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P12883
Hazard Information	This Myosin 7 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human breast carcin...