

Product Datasheet

Anti-Paxillin Pxn Antibody (Monoclonal, PXC-10) (orb18256)

Description	Anti-Paxillin Pxn Antibody (Monoclonal, PXC-10). Tested in ICC, WB applications. This antibody reacts with Chicken, Human, Mouse, Rat.
Species/Host	Mouse
Reactivity	Bovine, Gallus, Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	C-terminal part of recombinant chicken paxillin (amino acids 305-559).
Form/Appearance	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Immunocytochemistry , 1µg/ml, Human, mouse, rat, chicken, -Western blot, 1-2µg/ml, Human, mouse, rat, chicken. Add 1ml of PBS buffer will yield a concentration of 100ug/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone Number	PXC-10
Antibody Type	Primary Antibody
MW	65 kDa

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Uniprot ID**Q66H76****Expiration Date**

12 months from date of receipt.



Western blot analysis of 293t cell lysate (lane1), hela cell lysate (lane2), mcf-7 cell lysate (lane3), m231 cell lysate (lane4), jurkat cell lysate (lane5) using Paxillin antibody

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)