



Product Datasheet

Retinoblastoma Antibody / Rb1 (orb1825602)

Description	Recognizes a 105kDa phosphoprotein, identified as retinoblastoma (Rb) gene product. Its epitope is I...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human full-length protein was used as the immunogen for the recombinant Retinoblastoma antibody.
Storage	Aliquot the recombinant Retinoblastoma antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P06400
Hazard Information	This recombinant Retinoblastoma antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder carci...



SDS-PAGE analysis of purified, BSA-free ...