



Product Datasheet

Periostin Antibody (orb1825651)


Description	Periostin (PN), also designated osteoblast-specific factor 2 (OSF-2), is a disulfide linked protein ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 193-326) from the human protein was used as the immunogen for the recombinant POSTN antibody.
Storage	Aliquot the recombinant POSTN antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	Q15063
Hazard Information	This recombinant POSTN antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate tiss...



IHC staining of FFPE human lung tissue w...



SDS-PAGE analysis of purified, BSA-free ...