



Product Datasheet

Decorin Antibody (orb1825659)

Description	Decorin is a small leucine-rich proteoglycan (SLRP) family member that consists of a glycosaminoglyc...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant human Decorin protein fragment (within amino acids 212-336) was used as the immunogen for the Decorin antibody.
Storage	Aliquot the Decorin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P07585
Hazard Information	This Decorin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate carc...



SDS-PAGE analysis of purified, BSA-free ...