


Product Datasheet

Interferon gamma Antibody / IFNG (orb1825669)

Description	Recognizes a protein of 20-25kDa, identified as human interferon. This mAb is specific to human IFN- ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF
Immunogen	Purified recombinant human IFN-gamma was used as the immunogen for the Interferon gamma antibody.
Storage	Aliquot the Interferon gamma antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P01579
Hazard Information	This Interferon gamma antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml
Expiration Date	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...