



Product Datasheet

STAT6 Antibody (orb1825676)

Description	STAT6 is a transcription factor in the Jak/STAT signal transduction pathway responsible for mediatin...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 647-847) from the human protein was used as the immunogen for the STAT6 antibody.
Storage	Aliquot the STAT6 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P42226
Hazard Information	This STAT6 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human cerebellum ti...



IHC staining of FFPE human liposarcoma t...