



Product Datasheet

CD63 Antibody / LAMP-3 (orb1825686)

| | |
|----------------------------|---|
| Description | This mAb recognizes protein of 26kDa-60kDa, which is identified as CD63. Its epitope is different fr... |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P |
| Immunogen | Recombinant full-length human protein was used as the immunogen for the recombinant LAMP-3 antibody. Its epitope is different from that of mAb LAMP3/529. |
| Storage | Aliquot the recombinant LAMP-3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Formula | 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide |
| Isotype | Rabbit IgG, kappa |
| Clonality | Recombinant |
| Uniprot ID | P08962 |
| Hazard Information | This recombinant LAMP-3 antibody is available for research use only. |
| Dilution Range | Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT |
| Expiration Date | 12 months from date of receipt. |



IHC staining of FFPE human tonsil tissue...



SDS-PAGE analysis of purified, BSA-free ...