



Product Datasheet

CDK4 Antibody (orb1825694)


Description	Cyclin-dependent kinase-4 (CDK4) is a protein-serine kinase involved in the cell cycle. It is essent...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment corresponding to the N-terminal of human CDK4 protein was used as the immunogen for the recombinant CDK4 antibody.
Storage	Aliquot the recombinant CDK4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P11802
Hazard Information	This recombinant CDK4 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder cancer...



IHC staining of FFPE human serous ovaria...



IHC staining of FFPE human mammary cancer...