



Product Datasheet

Cytokeratin 5/6 Antibody (orb1825708)

Description	Cytokeratin 5 is a basic protein that is closely related to Cytokeratin 6. They share similar tissue...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Human Cytokeratin 5- and Cytokeratin 6-homologous recombinant protein fragments were used as the immunogen for the Cytokeratin 5/6 antibody.
Storage	Aliquot the Cytokeratin 5/6 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P02538, P13647
Hazard Information	This Cytokeratin 5/6 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate carc...



Negative control: IHC testing of FFPE hu...