


Product Datasheet

Basic Cytokeratin Antibody (HMW / Type II) (orb1825733)

Description	This mAb recognizes basic (Type II or HMW) cytokeratins, which include 67kDa (CK1); 64kDa (CK3); 59k...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human KRT76 protein was used as the immunogen for the recombinant Basic Cytokeratin antibody.
Storage	Aliquot the recombinant Basic Cytokeratin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2b, kappa
Clonality	Recombinant
Uniprot ID	Q01546
Hazard Information	This recombinant Basic Cytokeratin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human skin tissue w...