


## Product Datasheet

# Basic Cytokeratin Antibody (HMW / Type II) (orb1825735)

<b>Description</b>	This mAb recognizes basic (Type II or HMW) cytokeratins, which include 67kDa (CK1); 64kDa (CK3); 59k...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	Recombinant full-length human KRT76 protein was used as the immunogen for the recombinant Basic Cytokeratin antibody.
<b>Storage</b>	Aliquot the recombinant Basic Cytokeratin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
<b>Isotype</b>	Mouse IgG2b, kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>Q01546</b>
<b>Hazard Information</b>	This recombinant Basic Cytokeratin antibody is available for research use only.
<b>Dilution Range</b>	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
<b>Expiration Date</b>	12 months from date of receipt.



IHC staining of FFPE human skin tissue w...