



Product Datasheet

KRT7 Antibody / Cytokeratin 7 (orb1825737)

Description	It recognizes an intermediate filament protein (IFP) of 55kDa, which is identified as cytokeratin 7....
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human KRT7 protein was used as the immunogen for the recombinant KRT7 antibody.
Storage	Aliquot the recombinant KRT7 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Recombinant
Uniprot ID	P08729
Hazard Information	This recombinant KRT7 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate carc...



SDS-PAGE analysis of purified, BSA-free ...