



Product Datasheet

Cytokeratin 5 Antibody / KRT5 (orb1825744)


Description	This mAb recognizes a protein of 58kDa, which is identified as Cytokeratin 5 (KRT5). This type II cy...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 390-590) from the human protein was used as the immunogen for the recombinant KRT5 antibody.
Storage	Aliquot the recombinant KRT5 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Recombinant
Uniprot ID	P13647
Hazard Information	This recombinant KRT5 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate tissue...



IHC staining of FFPE human tonsil tissue...



SDS-PAGE analysis of purified, BSA-free ...