

Product Datasheet

FGF-23 Antibody / Fibroblast Growth Factor 23 (orb1825774)

Description

Fibroblast growth factor-1 (FGF-1), also designated

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

Applications

IHC-P

Immunogen

A recombinant partial protein sequence (within amino acids 1-251) from the human protein was used as the immunogen for the FGF23 antibody.

Storage

Aliquot the FGF23 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Uniprot ID

[Q9GZV9](#)

Hazard

Information


This FGF23 antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date


12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



SDS-PAGE analysis of purified, BSA-free ...



IHC staining of FFPE human kidney tissue...