



Product Datasheet

von Willebrand Factor Antibody / vWF (orb1825793)

Description	Human von Willebrand factor (VWF or factor VIII-related antigen) is a 309 kDa multimeric plasma glyco...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 1150-1250) from the human protein was used as the immunogen for the von Willebrand Factor antibody.
Storage	Aliquot the von Willebrand Factor antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P04275
Hazard Information	This von Willebrand Factor antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining
of FFPE
human
uterus
tissue...



SDS-PAGE
analysis of
purified,
BSA-free ...