



Product Datasheet

von Willebrand Factor Antibody / vWF (orb1825795)

Description	Human von Willebrand factor (VWF or factor VIII-related antigen) is a 309 kDa multimeric plasma glyco...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 1150-1250) from the human protein was used as the immunogen for the von Willebrand Factor antibody.
Storage	Aliquot the von Willebrand Factor antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P04275
Hazard Information	This von Willebrand Factor antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human uterus tissue...



SDS-PAGE analysis of purified, BSA-free ...