



## Product Datasheet

### CLEC4C Antibody / CD303 (orb1825817)

<b>Description</b>	Hyaluronic acid (HA) is a nonsulfated glycosaminoglycan that regulates cell adhesion and migration. ...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	A recombinant partial protein sequence (within amino acids 1-250) from the human protein was used as the immunogen for the CD303 antibody.
<b>Storage</b>	Aliquot the CD303 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	1 mg/ml in 1X PBS; BSA free, sodium azide free
<b>Isotype</b>	Mouse IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>Q8WTT0</b>
<b>Hazard Information</b>	This CD303 antibody is available for research use only.
<b>Dilution Range</b>	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
<b>Expiration Date</b>	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



IHC staining of FFPE human placental tis...