



Product Datasheet

FLI1 Antibody (orb1825927)


Description	Recognizes a protein of 51kDa, which is identified as FLI1. This protein, a member of the ETS family...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human FLI1 protein was used as the immunogen for the FLI1 antibody.
Storage	Aliquot the FLI1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q01543
Hazard Information	This FLI1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human Ewing's Sarco...



IHC staining of FFPE human Ewing's Sarco...



Western blot testing of human Daudi cell...