



Product Datasheet

FLI1 Antibody (orb1825946)

Description	Recognizes a protein of 51kDa, which is identified as FLI1. This protein, a member of the ETS family...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 200-400) from the human protein was used as the immunogen for the recombinant FLI1 antibody.
Storage	Aliquot the recombinant FLI1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q01543
Hazard Information	This recombinant FLI1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human lymph node ti...



SDS-PAGE analysis of purified, BSA-free ...