


Product Datasheet

SM22 alpha Antibody / Transgelin (orb1825959)

Description	This mAb recognizes a 22kDa protein, identified as Transgelin, also designated SM22-alpha. It may cr...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human Transgelin (TAGLN) protein was used as the immunogen for the recombinant Transgelin antibody.
Storage	Aliquot the recombinant Transgelin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Recombinant
Uniprot ID	Q01995
Hazard Information	This recombinant Transgelin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 0.5-1ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate carc...