



Product Datasheet

NCAM / CD56 Antibody (orb1825967)

Description	This mAb reacts with an extracellular domain (close to transmembrane) of CD56/NCAM. Three isoforms o...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 684-894) from the human protein was used as the immunogen for the recombinant NCAM antibody.
Storage	Aliquot the recombinant NCAM antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2b, kappa
Clonality	Recombinant
Uniprot ID	P13591 , P13592
Hazard Information	This recombinant NCAM antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon tissue ...



SDS-PAGE analysis of purified, BSA-free ...