



Product Datasheet

IL-18 Antibody (orb1825979)


| | |
|----------------------------|---|
| Description | Four structurally related IL-1 receptor ligands have been described. These include three agonists de... |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P |
| Immunogen | A recombinant partial protein (within amino acids 1-193) from the human protein was used as the immunogen for the IL-18 antibody. |
| Storage | Aliquot the IL-18 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Formula | 1 mg/ml in 1X PBS; BSA free, sodium azide free |
| Isotype | Mouse IgG2, kappa |
| Clonality | Monoclonal |
| Uniprot ID | Q14116 |
| Hazard Information | This IL-18 antibody is available for research use only. |
| Dilution Range | Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT |
| Expiration Date | 12 months from date of receipt. |



Analysis of a HuProt (TM) microarray con...



SDS-PAGE analysis of purified, BSA-free ...



IHC staining of FFPE human adrenal gland...