


Product Datasheet

NCAM / CD56 Antibody (orb1825996)

Description	This mAb reacts with an extracellular domain (close to transmembrane) of CD56/NCAM. Three isoforms o...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human NCAM1/CD56 protein was used as the immunogen for the NCAM antibody.
Storage	Aliquot the NCAM antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P13591, P13592
Hazard Information	This NCAM antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon tissue ...