


Product Datasheet

NCAM / CD56 Antibody (orb1826021)

| | |
|----------------------------|--|
| Description | This mAb reacts with an extracellular domain (close to transmembrane) of CD56/NCAM. Three isoforms o... |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Immunogen | A recombinant partial protein (within amino acids 400-650) from the human protein was used as the immunogen for the CD56 antibody. |
| Storage | Aliquot the CD56 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Formula | 1 mg/ml in 1X PBS; BSA free, sodium azide free |
| Isotype | Mouse IgG2a |
| Clonality | Monoclonal |
| Uniprot ID | P13591 , P13592 |
| Hazard Information | This CD56 antibody is available for research use only. |
| Dilution Range | Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT |
| Expiration Date | 12 months from date of receipt. |



Western blot testing of human brain tiss...



IHC staining of FFPE human colon tissue ...