



Product Datasheet

Progesterone Receptor Antibody (orb1826045)

Description	Progesterone is one of the central regulators of female reproduction. In breast development, progest...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 1-200) from the human protein was used as the immunogen for the Progesterone Receptor antibody.
Storage	Aliquot the Progesterone Receptor antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P06401
Hazard Information	This Progesterone Receptor antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human breast cancer...



SDS-PAGE analysis of purified, BSA-free ...